

J&Tec



Vetia

Anesthesia machine

Compact design for veterinary.
Low lifetime cost
anestentialator(option)

Vetia

Anesthesia machine

Technical specification

Safety

Meet the requirement of IEC60601 series, CE marking according to MDD93/42/EEC

Dimension

Dimension : 560mm(W) x 460mm(D) x 1250mm(H)

Operation Environment

Temperature : 0 ~ 40℃

Humidity : 30 ~ 85%

Gas supplying system

Range 1-gas single flowmeters

O₂ : 0 ~ 10LPM

Color coded safety system

Duble circulation CO2 absorber

Capacity Approximately 1kg x 1 canister

Patient gauge -20 ~ 100cmH₂O

Pressure relief valve 80cmH₂O

Pop-off valve

Selector valve (option)

Oxygen flush valve

self-returning push button

Flow 50LPM

Alarm whistle for oxygen supplying failure

Whistle sound at less than 2.5kgf/cm² of oxygen pressure

Vapor α temperature compensated vaporizer

Dial control set valve Halothane, Enflurane, Isoflurane : 0.2 ~ 6 Vol%

Sevoflurane : 0.2 ~ 8 Vol%

Temperature range 15 ~ 35 ℃

Flow range 0.5 ~ 15LPM

Mounting type Cage mount

Any Anesthetic Ventilator

(option)

General power supply : AC 85 ~ 264V, 50/60Hz

Battery : 12Vdc, 1.2AH

Option RS 232 serial data output

Display Rate, I:E ratio, Low/High pressure alarm limit, peak pressure

Option Volume monitor, Minute volume

Control Rate : 5 ~ 80 BPM

I:E ratio : 1:1 ~ 1:3

Tidal Volume : Max 1500ml

Low/High pressure alarm

Automatic and manual alarm setting for high/low alarm limit

Mute : 60sec

Stand-by

Audible and visual alarms Apnea, High/Low pressure alarm, Low battery



J&Tec Co., Ltd.

#139-2, Dongpae-ri, Gyoha-eup, Paju-si,
Gyeonggi-do, Korea
Tel : +82 31 957 7972 Fax : +86 31 942 1696
Email : JandTec@paran.com Website : www.JandTec.com

DISTRIBUTOR:

NP Intertrade Co.,Ltd.
121/110 Soi Saimai 21, Saimai Rd., Saimai, Bangkok 10220 Thailand
Tel : 662-991-3391 Fax : 662-991-3392 Mobile : 665-136-4688
E-mail : vetpook@rocketmail.com