


Innovative Usability for Enhanced Safety



Drive Gas Free
With electronically driven by high-performance turbine



8-Inch Capacitive Touch Screen
For intuitive operation and efficient interaction



Automatic Leak Test
Graphical guidance and one-click start



Inspiration Hold Function
Prevents respiratory motion interference with CT imaging and surgery

Maintain the plateau pressure

Prevent respiratory motion interference

Configurable start delay & hold duration



Veterinary Anesthesia Super Hub



* Configuration details subject to actual order specifications



LinkedIn Facebook Instagram

Q | Search Mindray Animal Medical

www.mindrayanimal.com
service@mindrayanimal.com

P/N:ENG-Veta 5Plus-210285X6P-20240802
©2024 Shenzhen Mindray Animal Medical Technology Co., Ltd. All rights reserved.



Veta 5Plus
Veterinary Anesthesia Machine

Redefining Veterinary Anesthesia Machines





Empowering Safety Elevating Confidence

Since its inception, the Veta series of veterinary anesthesia machines has consistently provided safer and more user-friendly anesthesia methods for clinical use.

The Veta 5Plus integrates the concept of active safety into veterinary anesthesia, innovatively developing the high-level ventilation mode—VS+*, which combines the experience of clinical experts with Veta's well-tested ventilation strategies to offer clinicians an intelligent anesthesia experience comparable to L3 assisted driving.

Innovative VS+ A New Era of Active Safety

Active Ventilation Adjustment
Corrects EtCO₂ by adjusting mechanical ventilation parameters

Instant Response
Synchronized adjustments guided by real-time data

Target EtCO₂ Control
Continuously tracks patient EtCO₂

VS+

Based on VS Mode
Developed from well-tested VS patterns

Gradual Adjustment
Rapidly adjusts when EtCO₂ exceeds the set range

Ensuring Active Safety

Multiple Safeguards
Adjusts within safety ranges for volume/pressure/frequency

Highest-Level Alarms
Activates if EtCO₂ is not corrected, alerting for immediate intervention



Boundary-Breaking Connectivity Anesthesia Super Hub

Interconnecting data between monitors and anesthesia machines, The Veta 5Plus directly reads the patient's EtCO₂ from monitor and applies the VS+ ventilation strategy

Classic VS Superior Safety in Ventilation

Breathing Synchronization
Intelligent trigger for reduced respiratory resistance

VS

Backup Ventilation
Effective prevention of hypoventilation with automatic backup ventilation

Smart Pre-set
Automatic adjustment of mechanical ventilation parameters based on patient weight



*Before using VS+, please thoroughly read the user manual and continuously monitor the depth of anesthesia and the patient's condition during use.